

SIZING GUIDE Two-Stage Breast Reconstruction

MENTOR® CPX™ Breast Tissue Expander / MENTOR® CPG™ Gel Breast Implant

STAGE 1: Selection of expander.*

STAGE 2: Selection of MENTOR® CPG™ Gel Breast Implant to use as a long-term implant after expansion is complete.

STAGE 1

CPX™2 Expander (Style 6100), CPX™4 Expander (Style 8100), and CPX™4 Expander with Suture Tabs (Style 9100) Breast Tissue Expander - Low Height						
Volume cc	Catalog # 6100	Catalog # 8100	Catalog # 9100	Width cm	Height cm	Proj. cm
250	354-6111	354-8121	354-9121	11,4	8,1	6,1
350	354-6112	354-8122	354-9122	12,7	9,4	6,5
450	354-6113	354-8123	354-9123	14,0	10,2	7,1
550	354-6114	354-8124	354-9124	15,0	10,9	7,4
650	354-6115	354-8125	354-9125	15,7	11,2	7,9
750	354-6116	354-8126	354-9126	16,5	11,9	8,1

CPX™2 Expander (Style 6200), CPX™4 Expander (Style 8200), and CPX™4 Expander with Suture Tabs (Style 9200) Breast Tissue Expander - Medium Height						
Volume cc	Catalog # 6200	Catalog # 8200	Catalog # 9200	Width cm	Height cm	Proj. cm
275	354-6211	354-8221	354-9221	10,7	9,3	6,2
350	354-6212	354-8222	354-9222	11,7	10,0	6,6
450	354-6213	354-8223	354-9223	12,7	10,8	7,0
550	354-6214	354-8224	354-9224	13,5	11,7	7,4
650	354-6215	354-8225	354-9225	14,6	12,6	7,6
800	354-6216	354-8226	354-9226	15,6	13,3	8,0

STAGE 1: Selection of expander.*

STAGE 1

CPX™2 Expander (Style 6200), CPX™4 Expander (Style 8200), and CPX™4 Expander with Suture Tabs (Style 9200) Breast Tissue Expander - Medium Height						
Volume cc	Catalog # 6200	Catalog # 8200	Catalog # 9200	Width cm	Height cm	Proj. cm
275	354-6211	354-8221	354-9221	10,7	9,3	6,2
350	354-6212	354-8222	354-9222	11,7	10,0	6,6
450	354-6213	354-8223	354-9223	12,7	10,8	7,0
550	354-6214	354-8224	354-9224	13,5	11,7	7,4
650	354-6215	354-8225	354-9225	14,6	12,6	7,6
800	354-6216	354-8226	354-9226	15,6	13,3	8,0

CPX™2 Expander (Style 6300), CPX™4 Expander (Style 8300), and CPX™4 Expander with Suture Tabs (Style 9300) Breast Tissue Expander - Tall Height						
Volume cc	Catalog # 6300	Catalog # 8300	Catalog # 9300	Width cm	Height cm	Proj. cm
250	354-6311	354-8321	354-9321	10,1	10,7	5,6
350	354-6312	354-8322	354-9322	11,3	11,8	6,0
450	354-6313	354-8323	354-9323	12,3	12,9	6,5
550	354-6314	354-8324	354-9324	13,2	13,8	6,9
650	354-6315	354-8325	354-9325	14,0	14,6	7,3
750	354-6316	354-8326	354-9326	14,6	15,3	7,6
850	354-6317	354-8327	354-9327	15,4	15,9	7,9

STAGE 2

CPG™ Implant 311 Low Height, Moderate Projection, Cohesive III™					CPG™ Implant 312 Low Height, Moderate Plus Projection, Cohesive III™					CPG™ Implant 313 Low Height, High Projection, Cohesive III™				
Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm
180	334-1101	11,0	9,8	3,7	225	334-1107	11,0	9,7	4,7	240	334-1106	11,0	9,9	4,9
210	334-1151	11,5	10,3	3,9	255	334-1157	11,5	10,2	4,9	270	334-1156	11,5	10,3	5,0
270	334-1251	12,5	11,2	4,2	330	334-1257	12,5	11,1	5,3	350	334-1256	12,5	11,2	5,4
300	334-1301	13,0	11,6	4,4	370	334-1307	13,0	11,5	5,5	395	334-1306	13,0	11,7	5,6
375	334-1401	14,0	12,5	4,7	465	334-1407	14,0	12,4	5,9	490	334-1406	14,0	12,6	6,1
460	334-1501	15,0	13,4	5,1	570	334-1507	15,0	13,3	6,3	605	334-1506	15,0	13,5	6,6
510	334-1551	15,5	13,8	5,2	690	334-1607	16,0	14,2	6,8					
560	334-1601	16,0	14,3	5,4										
615	334-1651	16,5	14,7	5,6										

CPG™ Implant 311 Low Height, Moderate Projection, Cohesive III™					CPG™ Implant 312 Low Height, Moderate Plus Projection, Cohesive III™					CPG™ Implant 313 Low Height, High Projection, Cohesive III™				
Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm
160	334-1051	10,5	9,4	3,6	195	334-1057	10,5	9,3	4,4	210	334-1056	10,5	9,5	4,8
180	334-1101	11,0	9,8	3,7	225	334-1107	11,0	9,7	4,7	240	334-1106	11,0	9,9	4,9
210	334-1151	11,5	10,3	3,9	255	334-1157	11,5	10,2	4,9	270	334-1156	11,5	10,3	5,0
235	334-1201	12,0	10,7	4,1	290	334-1207	12,0	10,6	5,1	310	334-1206	12,0	10,8	5,2
270	334-1251	12,5	11,2	4,2	330	334-1257	12,5	11,1	5,3	350	334-1256	12,5	11,2	5,4
300	334-1301	13,0	11,6	4,4	370	334-1307	13,0	11,5	5,5	395	334-1306	13,0	11,7	5,6
335	334-1351	13,5	12,0	4,6	415	334-1357	13,5	11,9	5,7	440	334-1356	13,5	12,2	5,8
415	334-1451	14,5	12,9	4,9	515	334-1457	14,5	12,8	6,1	545	334-1456	14,5	13,1	6,3
460	334-1501	15,0	13,4	5,1	570	334-1507	15,0	13,3	6,3	605	334-1506	15,0	13,5	6,6
510	334-1551	15,5	13,8	5,2	690	334-1607	16,0	14,2	6,8					
560	334-1601	16,0	14,3	5,4										

STAGE 2: Selection of MENTOR® CPG™ Gel Breast Implant to use as a long-term implant after expansion is complete.

STAGE 2

CPG™ Implant 321 Medium Height, Moderate Projection, Cohesive III™					CPG™ Implant 322 Medium Height, Moderate Plus Projection, Cohesive III™					CPG™ Implant 323 Medium Height, High Projection, Cohesive III™				
Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm
180	354-1058	10,5	9,9	3,8	225	334-1055	10,5	9,9	4,4	260	334-1052	10,5	9,9	5,3
215	354-1108	11,0	10,3	3,9	255	334-1105	11,0	10,3	4,7	300	334-1102	11,0	10,3	5,6
245	354-1158	11,5	10,8	4,0	295	334-1155	11,5	10,8	4,9	345	334-1152	11,5	10,8	5,8
280	354-1208	12,0	11,3	4,2	330	334-1205	12,0	11,3	5,1	390	334-1202	12,0	11,3	6,0
315	354-1258	12,5	11,8	4,4	375	334-1255	12,5	11,8	5,3	440	334-1252	12,5	11,8	6,2
355	354-1308	13,0	12,2	4,6	420	334-1305	13,0	12,2	5,5	495	334-1302	13,0	12,2	6,5
395	354-1358	13,5	12,7	4,7	475	334-1355	13,5	12,7	5,7	555	334-1352	13,5	12,7	6,7
480	354-1458	14,5	13,6	5,0	585	334-1455	14,5	13,6	6,1	685	334-1452	14,5	13,6	7,1
530	354-1508	15,0	14,1	5,2	650	334-1505	15,0	14,1	6,3					
640	354-1608	16,0	15,0	5,6										

CPG™ Implant 331 Tall Height, Moderate Projection, Cohesive III™					CPG™ Implant 332 Tall Height, Moderate Plus Projection, Cohesive III™					CPG™ Implant 333 Tall Height, High Projection, Cohesive III™				
Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm	Volume cc	Catalog #	Width cm	Height cm	Proj. cm
175	334-1003	10,0	10,2	3,4	205	334-1009	10,0	10,4	4,2	250	334-1004	10,0	10,4	5,1
200	334-1053	10,5	10,7	3,5										
230	334-1103	11,0	11,3	3,6	270	334-1109	11,0	11,5	4,7	330	334-1104	11,0	11,5	5,6
265	334-1153	11,5	11,8	3,7	305	334-1159	11,5	12,0	4,9	380	334-1154	11,5	12,0	5,8
300	334-1203	12,0	12,3	3,9	350	334-1209	12,0	12,5	5,1	430	334-1204	12,0	12,5	6,0
340	334-1253	12,5	12,8	4,0	395	334-1259	12,5	13,0	5,3	485	334-1254	12,5	13,0	6,2
380	334-1303	13,0	13,3	4,1	445	334-1309	13,0	13,5	5,5	545	334-1304	13,0	13,5	6,5
425	334-1353	13,5	13,8	4,3	495	334-1359	13,5	14,1	5,7	610	334-1354	13,5	14,1	6,7
475	334-1403	14,0	14,3	4,5	555	334-1409	14,0	14,6	5,9	680	334-1404	14,0	14,6	6,9
530	334-1453	14,5	14,8	4,6	615	334-1459	14,5	15,1	6,1	755	334-1454	14,5	15,1	7,1
585	334-1503	15,0	15,3	4,8	680	334-1509	15,0	15,6	6,3					
645	334-1553	15,5	15,9	5,0										

*PLEASE NOTE: Injection Dome is not MRI compatible.

For detailed indications, contraindications, warnings and precautions associated with the use of this device, please refer to the Product Insert Data Sheet.

Individual implant dimensions may vary slightly in products of this type. Not all units conform exactly to the dimensions noted above.